



**GLEN RESEARCH**

22825 DAVIS DRIVE  
STERLING, VIRGINIA  
20164

PHONE

703-437-6191

800-327-GLEN

FAX

703-435-9774

INTERNET

WWW.GLENRES.COM

## TAMRA NHS ESTER

### INTRODUCTION

Because rhodamine derivatives are generally not sufficiently stable to survive conventional deprotection, these can be attached to amino-modified oligonucleotides using post-synthesis labelling techniques. This product is the activated carboxylate, N-hydroxysuccinimide (NHS) ester, of tetramethyl rhodamine (TAMRA) dye in solution in dimethyl sulfoxide (DMSO) and is conjugated with an amino-modified oligonucleotide in sodium carbonate/bicarbonate buffer at pH 9. Although this technique is time consuming and places demands on the final purification to remove unconjugated dye, it is nevertheless routine and successful.

### USE OF TAMRA NHS ESTER

*Diluent:* ~0.17M Solution in Anhydrous DMSO

*Coupling:* Conjugate with an amino-modified oligo in sodium carbonate/bicarbonate buffer (pH=9).


*Storage:* Freezer storage, -10 to -30°C, dry

### LABELLING PROTOCOL

For a 0.2 μmole synthesis of an amino-modified oligo:

1. Dissolve oligo in 500 μL of conjugation buffer.
2. Add 6 μL of TAMRA/DMSO solution (~ 6 fold excess).
3. Vortex mixture and incubate at 37 °C in the dark for 1-2 hrs.
4. Separate oligo-TAMRA conjugate from salts and free TAMRA by size exclusion on a Glen Gel-Pak™ desalting column or equivalent.

**FIGURE 1: STRUCTURE OF TAMRA NHS ESTER**



Catalog Number: 50-5910-xx

### DESALTING USING A NAP-10 COLUMN

5. Equilibrate NAP column with approximately 10 mL of 50 mM TEAA buffer pH 7.
6. Load reaction mix on column and let flow into the column.
7. Add 0.5 mL TEAA buffer and let flow into column.
8. Elute oligo TAMRA conjugate with ≤ 1.5 mL TEAA buffer.
9. Collect the fast-running conjugate band and dry it in a vacuum concentrator.
10. Conjugate may be further purified by RP HPLC or PAGE to separate labelled from unlabelled oligonucleotides.

