

### EXTINCTION COEFFICIENTS AND FLUORESCENCE DATA

Calculate extinction coefficient of an oligo by either summing up the extinction coefficients of the individual bases times their number of occurrences. Or use a formula that takes into account nearest neighbor effects. An [algorithm for this calculation](#) can be found on the web. Just type in the sequence and the program will calculate the concentration of a  $\lambda$  A260/ml solution.

| Cat.No. | Nucleoside     | $\lambda_{\text{max-1}}$ | $E_{\text{max-1}}$ | $E_{\text{max-1}}$ | $\lambda_{\text{max-2}}$ | $E_{\text{max-2}}$ | E260             | E260       | Notes |
|---------|----------------|--------------------------|--------------------|--------------------|--------------------------|--------------------|------------------|------------|-------|
|         |                | (nm)                     | (ml/ $\mu$ mole)   | (L/mol.cm)         | (nm)                     | (ml/ $\mu$ mole)   | (ml/ $\mu$ mole) | (L/mol.cm) |       |
| 10-1001 | 7-deaza-dA     | 270                      | 11.3               |                    |                          |                    | 9.4              |            |       |
| 10-1003 | N6-Me-dA       | 266                      | 16.9               |                    |                          |                    | 15.2             |            |       |
| 10-1006 | Etheno-dA      | 295                      | 3.4                |                    | 274                      | 5.9                | 4.7              |            |       |
| 10-1007 | 8-Br-dA        | 266                      | 16.4               |                    |                          |                    | 14.8             |            |       |
| 10-1008 | 8-Oxo-dA       | 268                      | 12.2               |                    |                          |                    | 11.1             |            |       |
| 10-1014 | pdC            | 295                      | 7.7                |                    | 234                      | 14.7               | 5.1              |            |       |
| 10-1017 | Pyrrolo-dC     | 339                      | 2.36               |                    | 229                      | 17.5               | 2.41             |            |       |
| 10-1021 | 7-deaza-dG     | 259                      | 12.6               |                    |                          |                    | 12.6             |            |       |
| 10-1027 | 8-Br-dG        | 253                      | 12.1               |                    |                          |                    | 11.3             |            |       |
| 10-1028 | 8-oxo-dG       | 294                      | 5.2                |                    | 250                      | 6.7                | 5.9              |            |       |
| 10-1031 | 5'-OMe-dT      | 266                      | 9                  |                    |                          |                    | 8.3              |            |       |
| 10-1035 | Carboxy-dT     | 297                      | 16.1               |                    | 261                      | 14.7               | 14.7             |            |       |
| 10-1036 | 2-thio-dT      | 278                      | 17.5               |                    | 220                      | 14.8               | 10               |            |       |
| 10-1040 | dI (Inosine)   | 249                      | 12.5               |                    |                          |                    | 7.5              |            |       |
| 10-1041 | dNebularine    | 262                      | 7.1                |                    |                          |                    | 7                |            |       |
| 10-1043 | 3-Nitropyrrole | 283                      | 8.8                |                    |                          |                    | 7.7              |            |       |
| 10-1044 | 5-Nitroindole  | 328                      | 8.5                |                    | 265                      | 17                 | 16               |            |       |
| 10-1045 | 4-Methylindole | 265                      | 7.9                |                    |                          |                    | 7.2              |            |       |
| 10-1046 | 2-Aminopurine  | 303                      | 6.8                |                    | 243                      | 5.7                | 1                |            |       |
| 10-1047 | dP             | 294                      | 6.7                |                    | 231                      | 7.4                | 2.9              |            |       |
| 10-1048 | dK             | 279                      | 10.7               |                    |                          |                    | 7.7              |            |       |
| 10-1050 | dU             | 262                      | 10                 |                    |                          |                    | 10               |            |       |

|          |                     |          |      |        |     |      |      |                     |  |
|----------|---------------------|----------|------|--------|-----|------|------|---------------------|--|
| 10-1052  | 4-thio-dU           | 330      | 30.4 |        |     |      | 3.6  |                     |  |
| 10-1053  | 5-OH-dU             | 280      | 7.8  |        |     |      | 4.9  |                     |  |
| 10-1054  | pdU                 | 291      | 11.3 |        | 231 | 11.4 | 3.5  |                     |  |
| 10-1055  | d-pseudoU           | 262      | 7.7  |        |     |      | 7.6  |                     |  |
| 10-1056  | Fluorescein-dT      |          |      |        |     |      |      | 7.1<br>L/(mmol·cm)  |  |
| 10-1060  | 5-Me-dC             | 277      | 9    |        |     |      | 5.7  | 7.8<br>L/(mmol·cm)  |  |
| 10-1061  | 5-Me-dZ             | 314      | 4.8  |        | 218 | 8.6  | 1.8  |                     |  |
| 10-1062  | hmdC                |          |      |        |     |      |      | 8.7<br>L/(mmol·cm)  |  |
| 10-1063  | 5-OH-dC             | 292      | 6.3  |        | 220 | 13.3 | 3.4  |                     |  |
| 10-1065  | 5-Me-isodC          | 260      | 6.3  |        |     |      | 6.3  |                     |  |
| 10-1067  | 5-Me-isodC          | 260      | 6.3  |        |     |      | 6.3  |                     |  |
| 10-1076  | 7-deaza-dX          | 284      | 6.5  |        | 252 | 10.4 | 8.8  |                     |  |
| 10-1077  | iso-dG              | 292      | 11   |        |     |      | 4.6  |                     |  |
| 10-1078  | iso-dG              | 292      | 11   |        |     |      | 4.6  |                     |  |
| 10-1080  | 5-Br-dC             | 287      | 6    |        |     |      | 3.1  |                     |  |
| 10-1081  | 5-I-dC              | 293      | 5.7  |        |     |      | 3.3  |                     |  |
| 10-1085  | 2,6-diaminoPurine   | 278      | 10.2 |        | 255 | 9.3  | 8.5  |                     |  |
| 10-1090  | 5-Br-dU             | 278      | 9.7  |        |     |      | 5.1  |                     |  |
| 10-1091  | 5-I-dU              | 287      | 7.7  |        |     |      | 3.7  |                     |  |
| 10-1094  | Furano-dT           | See plot |      |        |     |      |      |                     |  |
| 10-1095  | 2,4-difluorotoluene | 266      | 2.3  |        |     |      | 1.8  |                     |  |
| 10-1097* | AP-dC               | 362      | 10.5 |        |     |      | 10.9 |                     |  |
| 10-1514  | Formyl-dC           |          |      |        |     |      |      | 11.3<br>L/(mmol·cm) |  |
| 10-1530  | dihydro-dT          | 210      | 6.3  |        |     |      | <0.1 |                     |  |
| 10-1539  | DBCO-dT             |          |      |        |     |      |      | 18,800              |  |
| 10-1550  | dihydro-dU          | 210      | 6.3  |        |     |      | <0.1 |                     |  |
| 10-1554  | 5-Ethynyl-dU        | 287      | 10.3 |        |     |      | 3.8  |                     |  |
| 10-1934  | 5-Formylindole      | 298      |      | 14,050 |     |      |      | 9,700               |  |
| 10-      | DBCO-TEG            |          |      |        |     |      |      | 8,000               |  |

|         |               |  |  |  |  |  |  |       |      |   |
|---------|---------------|--|--|--|--|--|--|-------|------|---|
| 1941    |               |  |  |  |  |  |  |       |      |   |
| 10-3055 | PseudoUridine |  |  |  |  |  |  | 8,100 | Ref: | <a href="http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3851715/">http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3851715/</a> |

Note: Biotin and Cholesterol have no absorbance at 260nm.

\*With an extinction coefficient of approximately 10,500 M<sup>-1</sup> and a quantum yield of fluorescence of 0.2, AP-dC is 2-3 times as bright as our popular Pyrrolo-dC analog. In addition, AP-dC exhibits a Stokes' shift greater than 100 nm. As with most fluorescent base analogs, it is substantially quenched upon forming a duplex. The quantum yield drops to 0.1 while gaining significant structure in the emission spectrum (Figure 4), making it an ideal probe of DNA structure.

#### FLUORESCENCE DATA

| Dye                       | E 260 nm                             | E λ max                              | Excitation max | Emission max | QY        | Notes  |
|---------------------------|--------------------------------------|--------------------------------------|----------------|--------------|-----------|--|
|                           | (M <sup>-1</sup> ·cm <sup>-1</sup> ) | (M <sup>-1</sup> ·cm <sup>-1</sup> ) | (nm)           | (nm)         |           |  |
| Acridine                  | 39,500                               | 9,120                                | 421            | 497          |           |  |
| 2-aminopurine             | 1,000                                | 3,600                                | 303            | 371          |           |  |
| 5'-CDPI <sub>3</sub> MGB™ | 37,900                               | 59,300                               | 340            |              |           |  |
| Cy3                       | 4,930                                | 136,000                              | 547            | 563          | 0.15      |  |
| Cy3.5                     | 24,000                               | 116,000                              | 591            | 604          | 0.15      |  |
| Cy5                       | 10,000                               | 250,000                              | 646            | 662          | 0.3       |  |
| Cy5.5                     | 21,500                               | 209,000                              | 688            | 707          | 0.3       |  |
| Dabcyl-dT                 | 29,100                               | 32,000                               | 476            |              |           |  |
| 5'-Dabcyl                 | 11,100                               | 32,000                               | 468            |              |           |  |
| Eclipse Quencher          | 6,600                                | 33,300                               | 530            | N/A          | 0         |  |
| Etheno-dA                 | 4,800                                | 5,800                                | 276            | 405          | 0.035     |  |
| Ferrocene-dT              | 14,200                               |                                      |                |              |           |  |
| 6-FAM                     | 20,900                               | 75,000                               | 495            | 521          | 0.9       |  |
| Fluorescein-dT            | 38,800                               | 75,000                               | 494            | 522          | 0.9       |  |
| HEX                       | 31,580                               | 96,000                               | 537            | 556          | 0.7       |  |
| Methylene Blue            | 10,300                               | 81,000                               | 665            |              |           |  |
| NBD                       | 3,700                                | 19,500                               | 485            | 535          | 0.1       |  |
| Psoralen                  | 16,500                               | 11,000                               | 301            |              |           |  |
| Pyrrolo-dC                | 4,000                                | 3,700                                | 345            | 470          | 0.07/0.02 | QY 0.07 single-stranded; 0.02 ds, deprotected in ammonia 55°C ON |
| Pyrene-dU                 | 18,500                               | 42,200                               | 402            |              |           |  |
| Redmond Red               | 12,100                               | 74,000 (pH 9.1)<br>52,300 (pH 7.1)   | 579            | 595          | 0.84      |  |
| TAMRA                     | 32,300                               | 89,000                               | 556            | 580          | 0.7       |  |
| TET                       | 16,255                               | 86,000                               | 519            | 539          | 0.9       |  |
| Yakima Yellow             | 23,700                               | 83,800                               | 530.5          | 549          | 0.96      |  |

#### PHYSICAL PROPERTIES OF BLACK HOLE QUENCHERS

| Quencher | λ <sub>max</sub> (nm) | E260 (L/mol.cm) | E <sub>max</sub> (L/mol.cm) |
|----------|-----------------------|-----------------|-----------------------------|
| BHQ-0    | 493                   | 7,700           | 34,000                      |
| BHQ-1    | 534                   | 8,000           | 34,000                      |
| BHQ-2    | 579                   | 8,000           | 38,000                      |
| BHQ-3    | 672                   | 13,000          | 42,700                      |

