


### 1. Product and Company Identification

<b>Product Code:</b>	<b>10-1891-xx</b>
<b>Product Name:</b>	Methacrylate C6 Phosphoramidite
<b>Product Use:</b>	For Research Use Only
<b>Manufacturer Information</b>	
<b>Company Name:</b>	Glen Research LLC 22825 Davis Drive Sterling, VA 20164 USA
<b>Phone Number:</b>	1-703-437-6191
<b>Web site address:</b>	www.glenresearch.com
<b>Email address:</b>	msds@glenres.com
<b>Information in case of emergency</b>	1-703-437-6191
<b>CAS Number:</b>	No data
<b>Synonyms:</b>	6-methacrylamido-hexyl-(2-cyanoethyl)-(N,N-diisopropyl)-phosphoramidite



### 2. Hazards Identification

#### GHS Hazard Statements

- H303: May be harmful if swallowed.
- H313: May be harmful in contact with skin.
- H320: Causes eye irritation.
- H333: May be harmful if inhaled.

#### GHS Precautionary Statements

- P103: Read label before use.
- P232: Protect from moisture.
- P233: Keep container tightly closed.
- P280: Wear protective gloves / protective clothing / eye protection / face protection.

#### GHS Response Statements

- P301 + 330: IF SWALLOWED: Rinse mouth with water.
- P302 + 350: IF ON SKIN: Gently wash with soap and water.
- P304 + 340: IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- P305 + 351 + 338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so - continue rinsing.
- P306 + 360: IF ON CLOTHING: Rinse contaminated clothing and skin immediately with plenty of water before removing clothes.
- P314: Get medical advice / attention if you feel unwell.
- P362: Take off contaminated clothing and wash before use.

#### GHS Storage and Disposal Phrases

- P401 Store in a cool, dry place
- P501: Dispose of contents / container in a safe way in accordance with all federal, state and local regulations.

**SAFETY DATA SHEET**  
**Methacrylate C6 Phosphoramidite**  
**(10-1891-xx)**

**Emergency Overview**

No known OSHA hazards. May cause eye, skin, and respiratory tract irritation.

**Route(s) of Entry:**

Inhalation? Yes      Skin? Yes      Eyes? Yes      Ingestion? Yes

**Potential Health Effects (Acute and Chronic)**

No data available.

**Recommended Exposure Limits**

No data available.

**LD 50 / LC 50**

No data available.

**Signs and Symptoms of Exposure**

No data available.

**Medical Conditions Generally Aggravated By Exposure**

No data available.

**3. Composition/Information on Ingredients**

**Essentially 100% Methacrylate C6 Phosphoramidite**

Product Name	CAS #	Concentration	M.W.	Formula
Methacrylate C6 Phosphoramidite	No data	100 %	385.48	C <sub>19</sub> H <sub>36</sub> N <sub>3</sub> O <sub>3</sub> P

**Synonyms:** 6-methacrylamido-hexyl-(2-cyanoethyl)-(N,N-diisopropyl)-phosphoramidite

The product contains no substance, which at their given concentration, are considered to be hazardous to health.

**4. First Aid Measures**

**Emergency and First Aid Procedures**

If inhaled, remove to fresh air.

If not breathing, give artificial respiration.

Get medical attention.

In case of skin contact, immediately wash skin with soap and plenty of water.

Remove contaminated clothing.

Wash clothing before reuse.

Get medical attention if symptoms occur after washing.

If swallowed, rinse mouth with water.

Get medical attention.

If in contact with eyes, rinse cautiously with water for several minutes.

Get medical attention if symptoms occur after washing.

Remove contact lenses if present and continue rinsing.

**SAFETY DATA SHEET**  
**Methacrylate C6 Phosphoramidite**  
**(10-1891-xx)**

**5. Fire Fighting Measures**

**Flash Point:** No data  
**Explosive Limits:** LEL: No data UEL: No data

**Fire Fighting Instructions**

A finely divided dust aerosol may be explosive, although the exact conditions of the hazards are unknown.  
Wear appropriate protective clothing to prevent contact with skin and eyes.  
Wear a self-contained breathing apparatus (SCBA) if necessary.

**Flammable Properties and Hazards**

No data available.

**Hazardous Combustion Products**

Oxides of Nitrogen, Carbon and Phosphorus

**Suitable Extinguishing Media**

Use a Class A Extinguisher (Dry chemical, carbon dioxide, water or foam).

**Unsuitable Extinguishing Media**

No data available.

**6. Accidental Release Measures**

**Steps To Be Taken In Case Material Is Released Or Spilled**

Sweep up bulk of material and finish by wiping with a damp towel. Do not let product enter the drain.  
Wear impermeable gloves, safety glasses and a lab coat when cleaning up the spill.  
A NIOSH dust mask is recommended when a possibility of dust inhalation exists. Local exhaust is usually sufficient.

**7. Handling and Storage**

**Precautions To Be Taken in Handling**

Handle using safe laboratory practices. Avoid all direct contact.

**Recommended Storage**

Freezer storage, -10 to -30°C, dry

**Precautions To Be Taken in Storing**

Keep container tightly closed. Store in well-ventilated place.

**Other Precautions**

Protect from sunlight.

**8. Exposure Controls/Personal Protection**

Product Name	CAS #	OSHA PEL	ACGIH TLV	Other Limits
Methacrylate C6 Phosphoramidite	No data	No data	No data	No data

**Respiratory Equipment (Specify Type)**

A NIOSH dust mask is recommended when a possibility of dust inhalation exists.

**Eye Protection**

Safety glasses

**Protective Gloves**

Impermeable gloves

**Other Protective Clothing**

Normal laboratory apparel

**Engineering Controls (Ventilation etc.)**

Local exhaust ventilation is usually sufficient

**Work/Hygienic/Maintenance Practices**

Wash hands after handling

**SAFETY DATA SHEET**  
**Methacrylate C6 Phosphoramidite**  
**(10-1891-xx)**

**9. Physical and Chemical Properties**

**Physical State:**  Gas  Liquid  Solid  
**Appearance:** Clear Viscous Syrup  
**Odor:** Odorless  
**Odor Threshold:** No data  
**pH:** No data  
**Melting Point:** No data  
**Freezing Point:** No data  
**Initial Boiling Point and Boiling Range:** No data  
**Flash Point:** No data  
**Specific Gravity:** ~1  
**Evaporation Rate (Water=1):** << 1  
**Flammability:** No data  
**Explosive Limits:** LEL: No data UEL: No data  
**Vapor Pressure (vs. Air or mm Hg):** No data  
**Vapor Density (vs. Air = 1):** No data  
**Relative Density:** No data  
**Solubility in Water:** No data  
**Partition Coefficient (n-octanol/water):** No data  
**Auto-ignition Temperature:** No data  
**Decomposition Temperature:** No data  
**Viscosity:** No data  
**Percent Volatile:** No data

**10. Stability and Reactivity**

**Stability:** Unstable  Stable   
**Conditions To Avoid - Instability**  
Stable under recommended storage conditions.  
**Incompatibility - Materials To Avoid**  
Strong oxidizing reagents  
**Hazardous Decomposition or Byproducts**  
Oxides of nitrogen, carbon and phosphorus may evolve when heated to decomposition.  
**Possibility of Hazardous Reactions:** Will occur  Will not occur   
**Conditions To Avoid - Hazardous Reactions**  
Stable material, hazardous polymerization will not occur.

**11. Toxicological Information**

To the best of our knowledge, the toxicological properties of this substance have not been investigated.  
This product should be handled with the usual care when dealing with chemicals.

Product Name	CAS #	NTP	IARC	OSHA
Methacrylate C6 Phosphoramidite	No data	No data	No data	No data

**12. Ecological Information**

No data available.

**13. Disposal Considerations**

**Waste Disposal Method**

For Product - Observe all federal, state and local regulations.  
For Contaminated Packaging - Dispose of in compliance with regulations.

**Methacrylate C6 Phosphoramidite****(10-1891-xx)**

Sweep up bulk of material and finish by wiping with a damp towel. Dispose of unused product.

**14. Transport Information****LAND TRANSPORT (49CFR)**

49CFR Shipping Name Not a dangerous substance as defined in the regulations.

**AIR TRANSPORT (ICAO/IATA)**

ICAO/IATA Shipping Name Not a dangerous substance as defined in the regulations.

**MARINE TRANSPORT (IMDG/IMO)**

IMDG/IMO Shipping Name Not a dangerous substance as defined in the regulations.

Marine Pollutant: No

**Additional Transport Information**

Not regulated as a hazardous material or dangerous goods.

**15. Regulatory Information****Regulatory Information**

SARA 302 Components: No chemicals in this product are subject to reporting requirements of SARA Title III, Section 302.

SARA 313 Components: This product does not contain any chemical components with known CAS numbers that exceed the threshold reporting levels established by SARA Title III, Section 313.

**16. Other Information**

For all other inquiries about this product contact Glen Research at 1-800-327-GLEN or 1-703-437-6191.

The information provided herein is based on sources believed to be reliable as of May 10, 2019 and pertains only to the material designated. Glen Research LLC makes no warranty or representation to its completeness, accuracy or currency. This material is intended for use by persons with the pertinent technical skills and at their discretion and risk. It is the responsibility of the user to determine the product's suitability for its intended use, the product's safe use and the product's proper disposal. Disposal of hazardous material may be subject to federal, state or local regulations.